SureFlex™ Laser Fiber

Reusable

Instructions for Use

SureFlex™ Laser Fiber <i>Reusable</i> Instructions for Use	page 1
Fibre laser SureFlex™ <i>réutilisable</i> Mode d'emploi	page 7
SureFlex™ Laserfaser <i>Wiederverwendbare</i> Gebrauchsanweisung	Seite 13
Fibra Laser SureFlex™ <i>Riutilizzabile</i> Istruzioni per l'uso	pagina 19
Fibra láser SureFlex™ <i>reutilizable</i> Instrucciones de uso	página 25
Fibra óptica de laser SureFlex™ <i>Reutilizável</i> Instruções de utilização	página 31
SureFlex™ laserfiber <i>Herbruikbaar</i> Gebruiksaanwijzing	pagina 37

Read All Instructions Thoroughly Before Use.

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Instructions for Use

DESCRIPTION

The SureFlex™ laser fiber is a fiber optic laser delivery device consisting of a patented, high energy termination (Black Hole™), a length of silica/silica fiber with a UV cured polymer secondary cladding and an ethylene tetrafluoro-ethylene copolymer (ETFE) jacket. The standard fiber output is a defect-free laser polish tip. SureFlex fibers are designed for use in a wide variety of surgical procedures as an integral part of a surgical laser system.

INDICATIONS

SureFlex laser fibers are intended for use in laser-based surgical applications, including, but not limited to endoscopic, laparoscopic and open surgical procedures involving vaporization, ablation, coagulation, hemostasis, excision, resection and incision of soft and cartilaginous tissue and surgical procedures involving vaporization, ablation and fragmentation of calculi (urinary and biliary). SureFlex is designed primarily for holmium (Ho:YAG) lasers, but may be used with any laser wavelength between 500nm and 2200nm that has been cleared for surgical use, e.g. Nd:YAG, KTP, Alexandrite, Argon, Ruby, Diode. SureFlex laser fibers may be reprocessed and reused up to five (5) times: detailed reprocessing instructions are provided within.

INSTRUCTIONS FOR USE

1. Prior to activation of the laser, review the laser Operator's Manual for instructions for proper set-up and operation of the laser system.
2. Remove the fiber from the sterile pouch using aseptic technique. The SureFlex™ fiber is provided sterile (EO) and may be reused if reprocessed and sterilized according to the instructions within. Do not use the product if the pouch seal is broken or if there is visual damage to the pouch.
3. Remove the fiber from the packaging starting from the distal end.

Do not start unwinding from the connector end.
Unwinding from the connector end will make fiber removal from the packaging difficult and may damage the fiber.
4. Make sure the laser is "OFF" or in "STANDBY" mode before inserting the fiber into the laser aperture (SMA port). Remove the protective cap from the connector and inspect the glass ferrule at the laser launch face for contamination. It is advisable to wipe the input face clean with a low lint wipe prior to use to remove accumulated organic contaminants. Attach

the connector to the laser and finger tighten the blue extension until a secure connection is achieved (generally indicated by release of lock-out mechanisms).

5. Handle the fiber with care: damage may occur if struck or bent sharply. Treat the output tip with particular care as it represents the most delicate, easily damaged portion of the assembly. If the fiber delivery system is damaged during use, it may be repaired (below) if sufficient length remains.
6. Activate the laser aiming beam and inspect the full length of the fiber for breaks. Breaks will appear as bright spots along the length of the fiber. When the tip is held at about 5 cm from a light colored surface, the aiming beam should be visible as a uniform, round spot. If the spot appears dim it may be useful to slightly loosen the connector nut and rotate the shaft of the connector within the laser port to optimize the aiming beam appearance, particularly in the case of small core fibers. If no breaks are visible and the spot is uniform, the fiber is ready for use.
7. Beginning treatment: in "free beam" applications, hold the probe about 5 cm from the distal tip during use. If using an endoscopic delivery system, insert the probe until the tip extends approximately 1 cm from the distal end of the endoscope or is in the field of view. Adjust the laser power to the level appropriate for the intended use and fiber size:
6-watts maximum for SureFlex™ 150, 12-watts maximum for SureFlex™ 200, 35-watts maximum for SureFlex™ 273, 100 watts maximum for SureFlex™ 365, 550, & 910.
8. Intraoperative cleaning - If the tip accumulates debris during use, turn the laser to "STANDBY" then carefully wipe the tip with a sterile gauze pad. Do not scrape the fiber with sharp instruments, as damage may result. Debris that is permitted to accumulate and "char" on the fiber tip will accelerate degradation in performance, reducing fiber lifetime and intraoperative reprocessing may be required.

INSTRUCTIONS FOR REPROCESSING AND REUSE

Fibers should be individually tracked to insure that they are not reused excessively. Steam sterilization exposure, particularly for fibers in tight coils, accelerates aging of the fiber. To insure appropriate safety in use, reuse should be limited to the approved maximum of five (5) cases. Before reuse, the SureFlex™ laser fiber must be cleaned and inspected for damage. The fiber may be stripped of buffer and cleaved for a fresh output face. Fibers must be sterilized before reuse.

CLEANING

1. Immerse the fiber for five (5) minutes (minimum) in an enzymatic detergent suitable for surgical instruments.
2. Use a soft gauze pad soaked in the detergent solution to wipe the fiber from end to end to remove all residue and debris.
3. Thoroughly rinse the fiber in warm tap water to remove the cleaning solution and any loosened residues.

4. Dry the fiber with a soft cloth, including the input face.

INSPECTION

1. With a magnifying eyepiece (loupe), microscope, endoscope or video system (minimum magnification 20X), inspect the fiber faces at both ends. While inspecting the fiber, slightly tilt the fiber face to view a variety of angles and reflected light from the fiber surface.
2. The fiber faces should be free of any defects such as chips, cracks, scratches, pits, blemishes or debris. Using a fiber with visible defects at the laser connector may damage the laser system. If defects are observed, do not reuse the fiber.
3. Degradation of the fiber face at the working tip is normal and will occur with use. If pitting, cracks, chips or debris are observed, the fiber should be stripped and cleaved before reuse.

CLEAVE AND STRIP THE DISTAL TIP OF THE FIBER

The recommended tools for cleaving and stripping fibers may be ordered from your authorized AMS distributor or directly from AMS. Check with an AMS sales representative for additional information. The following tools are recommended for stripping and cleaving the SureFlex™ laser fibers. See the accompanying instructions packaged with the Fiber Cleaver and the Fiber Stripper for information on use (including sterilization).

FIBER SIZE	TOOL P/N	GUIDE	COLOR
150µm Autoclavable	80-10001-003	0.0135"	Lt. Blue
200µm Autoclavable	80-10002-003	0.018"	White
273µm Autoclavable	80-10003-003	0.018"	Dk. Green
365µm Autoclavable	80-10005-003	0.031"	Yellow
550µm Autoclavable	80-10011-003	0.046"	Black
910µm Autoclavable	80-10008-003	0.063"	Dk. Brown

Cleaving device for all fibers: 10-0760

Ceramic SCISSORS ARE NOT RECOMMENDED

NOTE: Verify the stripping tool corresponds with the fiber size (tools for other manufacturers' fibers and other sizes of SureFlex™ fibers look similar and are often confused). Use the blade and buffer guard retainer color as a guide. Insure tools are sterile and maintained within the sterile field for intraoperative fiber cleaving and stripping.

STERILIZING THE FIBER

The fiber may be sterilized by EO or steam. The parameters shown below are offered **as a reference only**. Other validated sterilization methods and cycles may be used at the discretion of the safety and regulatory affairs officer. The fiber **MUST** be cleaned prior to sterilization. Optical fibers are glass and should be handled with some care. The individual user must validate sterility of reprocessed products prior to use.

EO STERILIZATION

1. Place LOOSELY coiled fiber in an EO compatible and permeable pouch and seal the pouch.
2. Precondition load (50 fibers maximum) at 80+% humidity for one hour at 130°F (54°C).
3. Sterilize in a 100% EO chamber with these conditions:
 - a) Gas Concentration: 700mg/l
 - b) Exposure Time: 1 hour
 - c) Temperature: 130°F (54°C)
 - d) Humidity: 80%
 - e) Aeration: 12 hours at 130°F (54°C)

STEAM STERILIZATION

1. Place LOOSELY coiled fiber in a steam sterilization pouch and seal the pouch.
2. Sterilize in a gravity displacement sterilizer steam sterilization chamber with these conditions:
 - a) Temperature: 270°F (132°C)
 - b) Exposure Time: 15 minutes
 - c) Drying Cycle: 15 minutes

PRE-VAC STEAM STERILIZATION

1. Place LOOSELY coiled fiber in a steam sterilization pouch and seal the pouch.
2. Sterilize in a forced air removal (pre-vac) steam sterilization chamber with these conditions:
 - a) # vacuum pulses: 3 pulses
 - b) Temperature: 275°F (135°C)
 - c) Exposure Time: 8 minutes
 - d) Drying Cycle: 25 minutes

NOTE ON FIBER PREPARATION FOR STERILIZATION: The tighter the fiber coil in sterilization (particularly in steam sterilization), the more fiber lifetime will be lost in the sterilization process. The cleaner the fiber input prior to sterilization, the longer it will remain transparent to the laser energy. How you handle the fiber has the most affect on how long the fiber lasts.

POTENTIAL COMPLICATIONS

Complications could include local and/or systemic infection, thermal damage to surrounding structures, local hematoma, dissection and perforation, tissue adhesion, distal tip detachment, and discomfort during and/or after energy application. In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of the tip. Properly reprocessed fibers, used within the accepted parameters of the field, should not lose tips. If tips detach, review reprocessing instructions.

CAUTIONS

- In soft tissues, increased laser exposure will result in a deeper and wider zone of necrosis.
- This product should only be used by physicians who have been trained in proper use and operation of surgical laser systems and are familiar with the physiological applications of laser energy.
- If tissue appears to be stuck to the lasing tip, put the laser in standby and carefully remove the tip from the target tissue. Remove the fiber from the scope. Use a soft, sterile cloth to clean the debris from the glass tip, taking care to avoid side pressure that may damage the tip.
- Do not bend fiber at sharp angles. If visible light (aiming beam) can be seen leaking from the fiber, fiber failure may result when therapeutic energy is applied as the fiber is deflected beyond the optical limits of total internal reflection.
- Care should be exercised with the glass tip to avoid severe impacts or side stresses that may fracture the tip.
- The time and energy required for treating a specific target area may be different when using the SureFlex™ laser fibers compared to using other fibers. A lower beam divergence, higher transmission efficiency and smaller laser spot size may alter treatment times.
- In soft tissue applications, during each lasing interval, it is helpful to move the distal tip back and forth or side to side 1 mm to 2 mm to prevent the tip from sticking to tissue.
- Immediately discontinue use if breaks or fractures appear in laser fiber. These breaks or fractures may allow undirected emission of laser energy, rendering the distal tip useless and potentially causing harm to surrounding tissues.

- **Always** insure all operating room personnel wear appropriate protective eyewear during the delivery of laser energy.
- Do **not** use the SureFlex™ laser fiber in the presence of flammable anesthetics or combustible materials.

HOW SUPPLIED AND STORAGE

- This SureFlex™ laser fiber is supplied sterile. Sterility is guaranteed only if the package is not opened or damaged and the product is used within the expiration date on the package label.
- Store in a cool, dry place.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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SureFlex™ and Black Hole™ are trademarks of AMS.
All rights reserved. Printed in USA.
American Medical Systems Canada Inc.
U.S.A.
American Medical Systems Europe B.V.
© 2014 American Medical Systems, Inc.